




# FIVALI HINGED KNEE BRACE

## Instructions for Use



### APPLICATION INFORMATION

1. Place the brace on the center of the kneecap, ensuring proper alignment with the patella (kneecap) ;
2. Wrap the bottom flap and then the top flap in front of the leg;
3. Secure the small inner band in a comfortable position;
4. Feed the top and bottom strap through the loop and then fix the strap with velcro strips;
5. Final application.

### CARE & WASHING INSTRUCTIONS

- Before the initial use of the knee brace, it is essential to wash it by hand and allow it to air dry naturally.
- Avoid the use of washing machines, as they can harm the structure and compression of the knitted fabric. Exercise caution to prevent the knee brace from snagging on other objects, as this can lead to damage.
- Ensure the product is not exposed to high-temperature environments, as this may cause scorching on the product's surface and result in a loss of elasticity. By following these guidelines, you can extend the lifespan and effectiveness of your knee brace.



### WARNINGS ⚠

- Neoprene supports should not be worn by individuals prone to dermatitis or allergic to neoprene products.
- If you experience pain, swelling, skin irritation or any other adverse reaction while using this product, consult your healthcare professional immediately.

### SINGLE USER ONLY

Please read the donning instructions and warnings carefully prior to use.  
To ensure that the brace works properly, follow all instructions.